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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUALITATIVE DETERMINATION OF IgG ANTIBODIES TO TOXOPLASMA GONDII IN CANINE AND FELINE SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the qualitative determination of IgG antibodies to *Toxoplasma gondii* in canine and feline serum or plasma.

This kit is designed for determination of IgG antibodies to *Toxoplasma gondii* in canine and feline serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 45 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Toxoplasmosis is a widespread infection caused by the intracellular protozoan parasite *Toxoplasma gondii*. Pets (cats and dogs) may serve as intermediate hosts for *T. gondii* and transmit infection to humans. In most cases, toxoplasmosis is a mild or asymptomatic disease; however, in immunocompromised patients this disease may be very severe and even life-threatening. Another risk group is pregnant women in whom primary toxoplasmosis can be transfected to the fetus, causing abortion or severe malformations. Testing of cats and dogs for toxoplasmosis is used to detect carrier state and thus to reveal the risk for humans.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by the antigen. Antibodies from the specimen bind coated antigen on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards species specific Ig, labeled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that infectious agents are not present in the reagents of this kit. All samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	Toxoplasma IgG-CF EIA strips, 8x12 wells	1	pcs		until exp.date
2 CAL	polystyrene microwells coated with antigen Toxoplasma gondii	1	pcs	blue	2 months
3 CONTROL – CONTROL +	feine or canine IgG antibodies to Toxoplasma gondii diluted in tris buffered preservative – 0,01% Bronidox L, 0,01% 2-Methyl-4-isothiazolin-3-one-hydrochloride; also contains blue dye dilution of preselected feline or canine serum, without IgG antibodies to Toxoplasma gondii with preservative – 0,01% Bronidox L, 0,01% 2-Methyl-4-isothiazolin-3-one-hydrochloride, colourless; dilution of preselected canine serum with high content of canine IgG antibodies to Toxoplasma Gondii; with preservative – 0,01% Bronidox L, 0,01% 2-Methyl-4-isothiazolin-3-one-hydrochloride; also contains red dye	2	pcs	colourless and red	2 months
4 CONJ HRP	aqueous solution of murine monoclonal antibodies to IgG coupled with horseradish peroxidase diluted on phosphate buffered solution with casein from bovine milk and detergent (Tween-20), contains 0,1% phenol as preservative and red dye	1	pcs	red	until exp.date
5 DIL SPE	phosphate buffered saline with casein from bovine milk and detergent (Tween-20), contains 0,1% phenol as preservative and blue dye	1	pcs	blue	until exp.date
6 SUBS TMB	ready-to-use single-component tetramethylbenzidine (TMB) solution.	1	pcs	colourless	until exp.date
7 BUF WASH 21X	aqueous solution of sodium chloride and detergent (Tween 20), contains proClin300 as a preservative	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 1 month at +2...+8°C or 5 days at RT
8 STOP	5,0% vol/vol solution of sulphuric acid	1	pcs	colourless	until exp.date
9 N003	Plate sealing tape	2	pcs		N/A
10 K101CFI	Instruction Toxoplasma IgG-CF EIA	1	pcs		N/A
11 K101CFQ	QC data sheet Toxoplasma IgG-CF EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl are useful but not essential;
- Calibrated micropipettes with variable volume, volume range 25–250 µl;
- Dry thermostat for 37 °C ±0,1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma. Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 fold dilution in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 6 wells for the calibrator CAL and control samples CONTROL -, CONTROL + and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute all samples using buffer DIL SPE (EIA buffer) 21 fold. (e.g. 25 µl of sample + 500 µl of diluent). After adding of samples (without heparin or EDTA), buffer changes colour from blue to red. Do not dilute control samples and calibrator.
3	Pipet 100 µl of calibrator, control samples CONTROL -, CONTROL + and diluted unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at 37 °C .
5	Prepare washing solution by 21 fold dilution of washing solution concentrate BUF WASH 21X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at 37 °C .
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells
10	Incubate 10–20 minutes at +18...+25 °C
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air
14	Apply coefficient method for data reduction. 1. Calculate mean OD for Calibrator. 2. Calculate a cut-off value by the following formulas: - for dogs: cut-off = OD(Calibrator) * Q value (please, see the QC insert) - for cats: cut-off = OD(Calibrator) * Q value *0.7 3. Divide mean OD of each sample by the relevant cut-off. $K = \text{OD}(\text{sample}) / \text{Cut-off}$ 4. If the K value is greater than 1.1, the result is POSITIVE. If the K value is less than 0.9, the result is NEGATIVE. If the K value is between 0.9 and 1.1, the result is EQUIVOCAL.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	Sample into the well, μl	Calculation factor
feline or canine serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	25 μl of sample + 500 μl of diluent	100	1

8. QUALITY CONTROL

It is recommended to use external control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a veterinarian to elaborate therapeutic measures. Each laboratory should establish its own normal range for Toxoplasma IgG-CF. Based on data obtained by XEMA, the following normal range is recommended (see below).

K-value more than 4.9 indicates high level response that may be associated with recent infection or vaccination.

Sex, age	Units, K	
	Lower limit	Upper limit
Seronegative	<0.1	0.9
Seropositive	1.1	4.9